



Take now **ALL THE TIME** you need for vascular suture, even in alarm condition of the patient

*Prenez dorénavant **TOUT LE TEMPS** nécessaire à la suture vasculaire, même en situation critique*

Dr Franceschi's HEMOSTATOR


The ideal Haemostasis Safety tool:

- > Immediate control of haemorrhage
- > Dramatic drop down of blood loss
- > Drying of operative field
- > Usable in every operative field even on non clamping zones (lumbar, pelvic vessels, sclerous tissues ...)

L'outil idéal d'hémostase de sauvetage :

- > Contrôle immédiat de l'hémorragie
- > Réduction majeure des pertes sanguines
- > Assèchement du champ opératoire
- > Utilisable quelque soit le site opératoire même sur des vaisseaux non clampables (lombaires, pelviens, ou au sein de tissus scléreux...)





The facts:

- > Emergency vascular suture is a very risky and stressful situation
- > Veins are fragile, less visible, and more difficult to dissect and clamp than arteries
- > Simply using the finger or dressing gauze that stop the bleeding hide the vessel
- > Pressure required to compress an aorta : only 1,5 kg

The Specifications:

- > 3 diameters for ideal fitting to the vessel size
- > The ring is lightly opened to allow the clearing of suture thread
- > The handle is oriented for a perfect control of the ring position

Les faits :

- > Une plaie vasculaire accidentelle provoque une situation stressante à haut risque
- > Les veines sont plus fragiles, moins visibles et plus difficiles à disséquer et clamer que les artères
- > L'utilisation du doigt ou d'une compresse pour arrêter le saignement cache le vaisseau
- > La pression nécessaire pour comprimer efficacement une aorte est seulement de 1,5 kg

Les specifications :

- > 3 diamètres pour s'adapter à la taille du vaisseau
- > Anneau discrètement ouvert pour permettre le dégagement du fil de suture
- > Poignée orientée pour un contrôle parfait de la position de l'anneau



Réf. DC53200 - 10



Réf. DC53200 - 20



Réf. DC53200 - 30